

ring User Guide - System Specifications & Flow

A Magic Window for Human Connection — Engineered by a team of 11 to 16-year-old girls

### [ Notice for Judges & Testers ]


To access the demo environments (Patient, Nurse, and Family portals), please use ID: 1 and Password: 1 on all login screens.

## 1For Patients: Zero-UI Technical Specifications

Even if physical movement is restricted, simply resting your gaze on the monitor ensures your thoughts are communicated. The "Zero-UI" eliminates complex physical interactions, utilizing the most natural human movement—the gaze—to drive the system.

### Hardware & Technical Specs (Patient Side)

- **Recommended Device:** Wall-mounted monitors or tablet devices (10 inches or larger).
- **Recognition Distance:** Optimal performance between 50cm to 80cm from the device camera.
- **Calibration:** AI maps pupil movement by having the patient look at corner points upon startup.
- **Dwell Time (Selection):** To prevent accidental inputs, an icon must be continuously stared at for 2.0 seconds to register as a "click."



## 2For Nurses: AI Triage & Social Experiment Data

---

Based on the reality that only 1.5% of calls are true emergencies, our AI triage mitigates alarm fatigue. Furthermore, "ring" serves as a platform for a large-scale social experiment beyond immediate care.

### Data & Technical Specs (Nurse & Admin Side)

- **Infrastructure:** Real-time data synchronization and triage powered by Firebase Cloud Functions.
- **Alert Categorization:** Instantly separates requests into "Emergency" (Red/Loud Alert) and "Non-Emergency/Chat" (Blue/Silent Notification).
- **Data Logging:** Records request types, timestamps, and response lead times. All data is completely anonymized for privacy.
- **Ultimate Goal:** Presenting this big data to the government to drive legal reforms for medical alert systems.



## 3For Families: Remote Care via Video Push

---

Providing the comfort of being in the same room, even 500km away. The family dashboard allows tracking of daily activities and features a crucial "Video Push" system to send warmth directly to the patient.

### Network & Technical Specs (Family Side)

- **Activity Sync:** View timelines of the patient's requests and summaries of their conversations with the AI.
- **Communication Protocol:** Real-time push notifications utilizing WebSockets.
- **Video Push:** MP4 videos uploaded from a family member's smartphone appear as a notification icon on the patient's screen.
- **Auto-Playback:** Simply by gazing at the notification, the video starts playing automatically in full screen with unmuted audio.

